



Neutron Camera NDCam

with New Generation MCP Neutron Detector



APPLICATIONS

- Nuclear Technology
- Material Research
- Non Destructive Testing
- Medical Science
- Biology
- Archeology/Paleontology
- Security

FEATURES

For more information visit

- 90° Angular Camera Setup
- 1.4 MPix CCD Camera (or better) out of Neutron Beam
- ❖ Based on ND 25 / ND 40 New Generation MCP Neutron Detector
- Active Area 25mm / 40 mm mm / Larger Size (coming soon)
- Excellent Cold and Thermal Neutron **Detection Efficiency**
- ❖ Spatial Resolution 50 µm (dependent upon neutron beam dispersion)

Version 1.1 / June 2014. Data may change without notice.