

T-COR 8™

The T-COR 8™ is the ideal point of care PCR instrument, offering all the features of a large powerful real-time PCR machine in a small package.

- 8 independently programmable amplification wells.
- 26.4 cm (10.4") touch screen with virtual keyboard.
- 4 hours of thermocycling battery life.
- Weighs less than 4.5 kg (9.9 lbs).
- Integrated barcode scanner.
- Custom protocol creator.
- Auto analysis with Smart CT™ software (calls positive/negative).
- Multiplex capable (FAM, CY3 and CY5).
- Run multiple protocols simultaneously.
- View data in real-time, chart or table format.
- View data from individual wells and channels.



- Data stored locally via SD card.
- Networkable via LAN for remote access and data download.
- Download data via CSV file through USB port.