


Simplicity

α/β Counter Control Software





Protean Instrument




Features

- One Software for All Protean Instrument Models
- Connects to up to 30 Separate Instruments
- Only Requires USB Connection
- Windows 7/Windows 10 Compatible
- SQL Express Database
- Alpha, Beta, and Guard* Plateau Options
- Quality Control Capabilities
- Simplified Sample Data Parameter Entry
- Data Exportable into Excel and CSV Formats
- Versatile Routine Programming Options
- Mass Curve Attenuation Capabilities
- Easy Sample Batch Generation
- Custom Report Editor
- Multi-level User Security Access
- ISO 11929 MDA Calculation Capable



* Guard Plateau not available on every model.



10744 Dutchtown Road, Knoxville, TN 37932 / <http://www.proteaninstrument.com>

Tel: 865-392-4600 / Fax: 865-392-4600 / E-mail: info@proteaninstrument.com